Laura Kangas, Milla Tynkkynen

& Mikko Korkiakoski BioBricks for measuring carbon starvation

Table of contents

Carbon starvation & how to measure it

System overview

Selected parts & details

Carbon starvation

 The cells undergo carbon starvation once they are running low on the carbon source (glucose)

- Insufficient carbon supply to metabolism causes a transition from the exponential growth phase to the stationary phase, reducing the reaction yield
- Carbon starvation in E. coli cells can lead to shrinkage of the bacterial cytoplasm

Measuring carbon starvation

- Low glucose concentrations in bioreactors lead to carbon starvation
- Low glucose concentration leads to the activation of adenylate cyclase, leading to formation of cAMP
- The PcstA promoter is activated by the binding of cAMP to the cAMP receptor protein (CRP) upstream of the transcription site, leading to transcription of the gene
- Thus, when glucose is not present, GFP will be produced, and glucose acts as a sort of repressor on transcription



System overview – Biological NOT gate



SELECTED PARTS

В В

BBa_K118011: PcstA, a glucose-repressible promoter

ğ

X

BBa_B0015: double terminator containing B0010-B0012



Chassis: *E. coli* strain BL21(DE3)

BBa_E0040: coding GFP



BBa_K1362091: plasmid pSB1A30, a high copy BioBrick cloning/expression backbone carrying Amp resistance



Measurements: fluorescent microscopy (visualization of GFP location and expression) and flow cytometry (fluorescent intensity)

DETAILS OF THE SELECTED PARTS

- BBa_E0040: coding GFP, gene of interest
 - Green fluorescent protein gene derived from jellyfish Aequeora Victoria
 - Excitation wavelength 545 nm (light absorbation)
 - Emission wavelength 475 nm (light emission)
 - GFP chosen to easily detect the expression
- BBa_B0015: double terminator containing B0010-B0012
 - B0010: transcriptional terminator derived from E.
 coli, 64 bp stem-loop
 - B0012: transcriptional terminator from coliphage T7, for E. coli RNA polymerase, promoter in the reverse direction
 - very commonly used terminator

DETAILS OF THE SELECTED PARTS

- BBa_K1362091: plasmid Psb1a30
 - a high copy BioBrick cloning/expression backbone
 - contains an Ampicillin resistance gene

MEASUREMENTS

- fluorescent microscopy
 - visualization of GFP location and expression
- flow cytometry

- laser used to detect and analyze cells
- measuring the intensity of fluorescence

Thank you!



Shi, H., Westfall, C. S., Kao, J., Odermatt, P. D., Anderson, S. E., Cesar, S., Sievert, M., Moore, J., Gonzalez, C. G., Zhang, L., Elias, J. E., Chang, F., Huang, K. C., & Levin, P. A. (2021). Starvation induces shrinkage of the bacterial cytoplasm. *Proceedings of the National Academy of Sciences of the United States of America*, *118*(24), e2104686118. <u>https://doi.org/10.1073/pnas.2104686118</u>)

Part:BBa K4115003 - parts.igem.org. Igem.org. Retrieved March 20, 2023, from http://parts.igem.org/wiki/index.php/Part:BBa_K4115003

Wei T, Dai H. Quantification of GFP signals by fluorescent microscopy and flow cytometry. Methods Mol Biol. 2014;1163:23-31. doi: 10.1007/978-1-4939-0799-1_3. PMID: 24841297.